



Configuration

Frame with height and backrest adjustment

Castors with brakes 100 mm

Seat W44 x D45 cm with an opening front or rear from Ø12 cm

Backrest with canvas

Armrests with armpads

Legrests with height adjustable footrests

Insertion rails

Foldable bucket

Options

Suitcase for travelling

Upholstered footrests

Calf strap

Alternative options

Traveller base frame only with leg supports, compatible with Life, Life Special, Junior, Junior Special, or Trip seat inserts

Application possibilities

- Shower use
- Toilet use (roll-over)
- Travelling (disassemblable)
- Independent transfers
- Pressure ulcer prevention
- Incontinence prevention
- Reduced torso stability
- Hip angle restriction
- Amputation
- Varying leg length

Width wheelchair incl. wheels 54 cm

Seat width 44 cm

Seat height floor – top of the seat 53-60 cm

Length wheelchair ca. 86 cm

Weight wheelchair ca. 18 kg

Loadable weight 120 kg

Recommended distance:
wall to toilet rim From 45 cm

Toiletheight 40-50 cm

Weight suitcase 5 kg

Measurements suitcase H81 x B55 x D35 cm

These specifications are only an indication and strongly depend on configuration and adjustment.

Scan for the pre-fit form

