



Configuration

Frame with backrest adjustment
100/125 mm castors with brakes
Custom seat with an opening left/right/front/back/middle or no opening
Backrest with canvas or upholstered and/or side support
Push handles
Armrests in different shapes and heights
Custom upholstered armpads
Height-adjustable, (comfort) legrests with folding footrests
Height-adjustable one-piece footrest
Upholstered footrests
Adjustable headrest
Pelvic belt
Abduction bobbin
Calf strap/calf protector
Insertion rails for bedpan/toiletbucket
Stainless steel bedpan/toiletbucket
Other individual customisation

Application possibilities

- Shower use
- Toilet use (roll-over)
- Hip width 48-54 cm
- Pressure ulcer prevention
- Incontinence prevention
- Reduced torso stability
- Hip angle restriction
- Scoliosis
- Small or large knee angle
- Short or varying leg length
- Small or large ankle angle (equinus foot)
- Amputation
- Differing body measurements
- Stationary use (without castors)
- Tripping propulsion
- Individual customisation

Width wheelchair incl. wheels	65 cm
Seat width	48-50-52-54 cm
Seat height floor - top of the seat	48-69 cm
Length wheelchair ca.	88 cm
Weight wheelchair ca.	17 kg
Loadable weight	150 kg
Recommended distance: wall to toilet rim	From 55 cm
Toiletheight	40-50 cm

These specifications are only an indication and strongly depend on configuration and adjustment.

Scan for the pre-fit form

